

BIOCLEAN™

CROC™

Antistatic Cleanroom Shoe

BCFW

About

Comfortable, lightweight and durable the BioClean Crocs™ are constructed from Croslite™ closed-cell resin. Antistatic and slip-resistant they also have anti-microbial properties to resist bacteria and odor and have full body shock absorbing benefits.

Specifications

MATERIAL:

Croslite™ closed-cell resin

MATERIAL:

Croslite™ closed-cell resin

COMPATIBILITY:

BSEN ISO 20347:2004, ASTM F 2413:2005.

COMPATIBILITY:

BSEN ISO 20347:2004, ASTM F 2413:2005.

CONSTRUCTION:

Croslite™ closed-cell resin

CONSTRUCTION:

Croslite™ closed-cell resin

COLOUR:

White

Features

- ESD tested to ASTM F2413:2005
- Anti-microbial properties to resist bacteria and odor
- Slip-resistant
- Ergonomic design proven to reduce fatigue
- Full body shock absorbing benefits
- Heel strap keeps shoe on foot, or rolled forward to create a slip-on

ANTISTATIC & SLIP-RESISTANT

ERGONOMIC DESIGN

ANTI-MICROBIAL PROPERTIES

SHOCK ABSORBING



Quality Standards

- Conforms to ASTM F1677 slip-resistance standards
- BS EN ISO 20347:2004 and ASTM F2413:2005 compatible
- Point to ground resistance 1×10^6 - 1×10^8 ohms (at 100 volts) as per ASTM F 2413:2005
- Typical Particle Count: <750 particles/cm² (Particle size: $>0.5\mu\text{m}$)

Shelf Life & Storage

Store in a dry, cool place ($<40^\circ\text{C}$) away from direct sunlight and fluorescent light.

TO ORDER

Re-order code	Size
BCFW-M2	EU 34
BCFW-M11	EU 44
BCFW-M12	EU 45

Packing

One pair per sealed inner PE bag, one inner bag per sealed outer PE bag; 12 outer bags per lined carton (12 pairs)

Note: To clean the Crocs™ use any concentration of soap/water solution, ensuring the temperature of the solution does not exceed 140°F. Cleaning may include using ESD Wipes, baby wipes or even a simple damp cloth soaked in water/soap/laundry detergent. Crocs™ can be wiped down with 70% IPA with no adverse effect.

Although the Crocs™ material is a durable material, it is susceptible to extreme heat, so the following should not be used in any cleaning regime;

Conventional dryers or dishwashers

Washing machines (with heat setting over 140°F)

Autoclave sterilization

CONTACT

Europe

(Headquarters)

Nitritex Ltd.
Minton Enterprise Park,
Oaks Drive,
Newmarket, Suffolk,
CB8 7YY,
United Kingdom
T: +44 (0) 1638 663338
F: +44 (0) 1638 668890
E: europa@bioclean.com

Other Regions

Asia:

malaysia@bioclean.com

Americas:

americas@bioclean.com



BioClean™ is a brand of Nitritex Ltd.

For more information, please visit www.bioclean.com

BCFW/1706/PDS9