

BIOCLEAN VECTOR

Non-sterile Vinyl Gloves
BVA-E

About

BioClean Vector Vinyl Gloves offer an economical alternative to nitrile whilst still providing tactility for precision work. The clear 300mm (12") long powder-free vinyl cleanroom gloves are ambidextrous and feature a beaded cuff for strength and stability on the arm.

Specifications

COMPATIBILITY:

ISO Class 5

LENGTH:

300mm (12")

MATERIAL:

Vinyl (PVC)

PROTEIN LEVEL:

PVC contains no natural latex proteins

SURFACE:

Smooth

SHAPE:

Ambidextrous

COLOUR:

Clear

NON-TEXTURED

LATEX-FREE

TACTILE

Features

- Non-textured
- Beaded cuff for strength
- Offer good tactility



Quality Standards

- Manufactured in a facility holding ISO 9001:2015 and (EU)2016/425 Module D certifications
- Conforms to Category 1 - Simple Design Personal Protective Equipment PPE Regulation (EU) 2016/425
- Processed in an NEBB certified ISO Class 5 environment

Physical Properties

CHARACTERISTICS	VALUE	TEST METHOD
Freedom From Holes	1.5 AQL	EN 374-2:2015
Glove Length +/-10mm (mm/inches)	290/11.4	EN 420:2003 + A1:2009
Single Wall Thickness Minimum (mils/mm)	Cuff: 2.36/0.06 Palm: 3.54/0.09 Finger Tip: 3.94/0.10	EN 455-2:2015
Force At Break During Shelf Life	≥3.6 N	EN 455-2:2015

PALM WIDTHS (mm) ± 5mm				
SIZE	Small	Medium	Large	XLarge
mm	85	95	105	115
TEST METHOD	EN 420:2003 + A1:2009			

Shelf Life & Storage

Five (5) years from date of manufacture.
Store in a dry, cool place (<40°C) away
from direct sunlight and fluorescent
light.

Glove Cleanliness Characteristics

PARTICLES				TYPICAL PARTICLE COUNT				TEST METHOD
$\geq 0.5\mu\text{m}(\text{counts}/\text{cm}^2)$				<3000				IEST-RP-CC005.4
TYPICAL EXTRACTABLE IONS								TEST METHOD
ANIONS	Fluoride	Chloride	Nitrite	Bromide	Nitrate	Phosphate	Sulphate	IEST-RP-CC005.4
$(\mu\text{g}/\text{cm}^2)$	0.003	0.07	0.001	0.002	0.03	0.019	0.12	
CATIONS	Lithium	Sodium	Ammonium	Potassium	Calcium	Magnesium	Zinc	
$(\mu\text{g}/\text{cm}^2)$	0.001	0.22	ND	0.007	0.035	0.002	ND	
ND = Not Detected, NT = Not Tested								

TO ORDER

Re-order code

BVA-E-S
BVA-E-M
BVA-E-L
BVA-E-XL

Size

Small
Medium
Large
XLarge

Packing

100 pieces per sealed inner PE bag; one inner PE bag per sealed outer PE bag; 10 outer bags per lined carton (1000 pieces)

For additional information visit us at www.bioclean.com, or contact us at;

Nitritex Ltd

Minton Enterprise Park
Oaks Drive
Newmarket, Suffolk
CB8 7YY
United Kingdom
Tel: +44 1638 663338

Ansell US

111 Wood Avenue South
Suite 210
Iselin, NJ 08830
USA
Tel: +1 800 800 0444

Nitritex (M) Sdn.Bhd

2 Jalan Jurunilai U1/20
Seksyen U1, Hicom Glenmarie Industrial Park
40150 Shah Alam, Selangor
Malaysia
Tel: +603 5569 3857/3859

BVA-E/1904/PDS24



Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.

*Warning: No glove provides complete protection against all chemicals. Users must test the gloves against the particular chemicals and environment in which they will be used.

**Please see product validation pack or contact Ansell customer service for specific data on use of gloves with cytotoxic drugs. Gloves used for protection against such drugs must be selected specifically for the type of chemicals used.

Ansell, ® and TM are trademarks owned by Ansell Limited or one of its affiliates. © 2012 – 2019 Ansell Limited. All Rights Reserved.