

# BIOCLEAN CHEMPREP Sterile Prep Mat S-BCPM

## About

The BioClean ChemPrep™ Prep Mat enhances workbench safety for technicians, nurses and pharmacists when handling or administering chemotherapy drugs. The Prep Mat is a highly-absorbent, three-layered preparation mat consisting of two layers of polyethylene with a central layer of air-laid paper. The upper layer allows even distribution of liquid, the middle layer gives absorption properties and the bottom layer provides a defence against leakages.

## Specifications

### MATERIAL:

Top layer - Polyolefin apertured film

Middle layer - Air-laid paper, Bottom

Layer - Polyethylene

### WEIGHT:

210gsm

### LENGTH:

56cm (22") or 40cm (16")

### COMPATIBILITY:

ISO Class 5 & EU GMP Grade A

### COLOUR:

Blue

### STERILIZATION:

Gamma irradiation

## Features

- Latex-free 3-layer construction
- Mesh upper layer for even distribution
- Highly absorbent middle layer
- Impermeable bottom layer

HIGHLY ABSORBENT

LOW-LINTING

3-LAYER CONSTRUCTION



## Quality Standards

---

- Manufactured in a facility holding ISO 9001:2015
- Processed in an NEBB certified ISO Class 5 environment

## Sterilization

---

Method: Gamma irradiation  
Minimum Dose: 25kGy  
Sterility Assurance Level:  $10^{-6}$

## Shelf Life & Storage

---

Three (3) years from date of manufacture. Store in a dry, cool place (<40°C) away from direct sunlight and fluorescent light.

## Physcial/General Properties

---

Property Measured	Method Used	Units of Measurement	Typical Value
Sorptive Capacity	IEST-RP-CC004.3	ml/m <sup>2</sup>	852
Sorptive Capacity	IEST-RP-CC004.3	ml/prep mat	83.3



## TO ORDER

### Re-order code

### Size

S-BCPM-4025B

40cm x 25cm

S-BCPM-5640B

56cm x 40cm

### Packing

S-BPCM-4025B:

One piece per sealed inner PE bag; one inner bag per second PE bag; 50 double bags per sealed outer PE bag; four outer bags per lined carton (200 pieces)

S-BPCM-5640B:

One piece per sealed inner PE bag; one inner bag per second PE bag; 25 double bags per sealed outer PE bag; four outer bags per lined carton (100 pieces)

For additional information visit us at [www.bioclean.com](http://www.bioclean.com), or contact us at;

#### Nitritex Ltd

Minton Enterprise Park  
Oaks Drive  
Newmarket, Suffolk  
CB8 7YY  
United Kingdom  
Tel: +44 1638 663338

#### Ansell US

111 Wood Avenue South  
Suite 210  
Iselin, NJ 08830  
USA  
Tel: +1 800 800 0444

#### Nitritex (M) Sdn.Bhd

2 Jalan Jurunilai U1/20  
Seksyen U1, Hicom Glenmarie Industrial Park  
40150 Shah Alam, Selangor  
Malaysia  
Tel: +603 5569 3857/3859

S-BCPM/1811/PDS20



Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of accessory for a specific application.

Ansell, ® and TM are trademarks owned by Ansell Limited or one of its affiliates. © 2012 – 2020 Ansell Limited. All Rights Reserved.